

## TECHNICAL DATA SHEET

### Healthy Hooves Hoof Slipper HH 107

Revision Date: 18/02/2020

#### **Product Description**

Healthy Hooves® Hoof Slipper HH 107 is an easy and rapid fitting product to keep hoof care products in place, on the specified area of the hoof. The Hoof Slipper is a bandage with a spongy filling, this needs to be placed between the claws of the hoof and fastened by Velcro.

#### **Benefits**

1. Easy and rapid fitting.
2. Sponge takes care of ongoing supply of the; Healthy Hooves Hoof Care Plus Paste & Hoof Care Plus Spray.
3. Prevents contamination via manure.
4. Easy to refill.

#### **Properties**

Form:	Solid
Specific Gravity:	N/A
Colour:	Black and White
Odour:	Odourless
Solubility:	N/A
pH (Neat):	N/A
pH 1% in solution H <sub>2</sub> O	N/A

#### **Method of Use**

The narrow part of the Hoof Slipper is placed between the claws of the hoof, then the broader part is placed against the bulb of the heel. Once applied, the cow can then walk normally and any contamination via manure will be prevented.

The Spot Hoof Liquid and/or Paste will be pushed to the bulb of the heel.

The Hoof Slipper can easily be refilled during application, this can be done daily in the milking parlour.

#### **Process of Use**

Healthy Hooves® Hoof Slipper HH 107 should be used as follows for the best results:

Can be used in conjunction with the Healthy Hooves Hoof Care Plus Paste & Hoof Care Plus Spray.

#### **Packaging**

Pack of 10.